

微量冷凍型離心機 / MICRO REFRIGERATED CENTRIFUGE

■ CN-1000R



- LCD數位面板，具溫度及轉速、離心力雙重切換設定及顯示功能。
- 溫度範圍：4°C ~ 40°C。
- 10組記憶系統，可預先設定常用參數。
- 微電腦控制系統，操作容易。
- 不平衡自動斷電安全裝置。
- 運轉停止時，自動蜂鳴告示。
- 高性能高效率控制系統。
- 過熱保護及過載安全裝置。
- 不銹鋼製內桶，防腐蝕、耐衝擊。
- 故障指示及保險絲安全裝置。
- 自動緩衝加速、減速功能。
- 自動氣流冷卻循環。
- 低噪音結構設計，運轉安靜。
- LCD, temperature, speed and RCF setting and display.
- Temperature range: 4°C ~ 40°C .
- 10 memory programs for preset parameters.
- Microprocessor control for easy operation.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- Overheat and overload safety protection.
- Stainless steel inside bowl.
- Breakdown indication and fuse safety device.
- Auto accelerating and decelerating functions.
- Auto-cooling air circulation.
- Auto current-restraining function.
- Low noise.

Model	Speed (rpm)	Rotor	Max. RCF (x g)	Max. Capacity (ml)	Timer (min.)	Weight (kg)	Size (cm) W x D x H
CN-1000R	14,000	REA-1518	16000	1.5 / 2.0 x 18	99:59	35	55 x 42 x 28

AC 110 / 220V , 50 / 60Hz

古典基本型離心機 / CLASSIC TYPE CENTRIFUGE

■ C-12B



6ml x 4pcs



10ml x 4pcs



- 適用於實驗室及臨床檢驗。
- 可調速、操作容易。
- 6ml/10ml之試管套可供選擇。
- Suitable for laboratory and clinical tests.
- Adjustable speed design.
- Can choose with 6ml or 10 ml adapter .

Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Weight (kg)
C-12B	0~3,000	920	10 x 4	3
		820	6 x 4	

Power: AC 110 / 220V , 50 / 60Hz

Size: W244 x D244 x H250 (mm)

微量型離心機 / MICRO TYPE CENTRIFUGE

■ MCD-2000



- 數位式轉速及時間顯示。
- 短暫離心功能。
- 電動安全門扣、運轉中自動關閉。
- 運轉停止時，自動蜂鳴告示。
- ABS塑鋼外殼並附觀測視窗。
- 自動氣流冷卻循環。
- LED speed and time display.
- The transient centrifugal functions.
- Auto cover lock, auto power off when open.
- Auto buzzer when finish.
- ABS plastics outer covering ; see-through window.
- Auto-cooling air circulation.



Model	Rotor	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)
MCD-2000	REA-1518	14,000	16,000	1.5 / 2.0 x 18	99 : 59	8

Power: AC 110 / 220V , 50 / 60Hz

Size: W290 x D320 x H245 (mm)