

古典基本型離心機 / CLASSIC TYPE CENTRIFUGE

■ C-12B



6ml x 4pcs



10ml x 4pcs



- 適用於實驗室及臨床檢驗。
- 可調速、操作容易。
- 6ml/10ml之試管套可供選擇。
- Suitable for laboratory and clinical tests.
- Adjustable speed design.
- Can choose with 6ml or 10 ml adapter .

Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Weight (kg)
C-12B	0~3,000	920	10 x 4	3
		820	6 x 4	

Power: AC 110 / 220V , 50 / 60Hz

Size: W244 x D244 x H250 (mm)

微量型離心機 / MICRO TYPE CENTRIFUGE

■ MCD-2000



- 數位式轉速及時間顯示。
- 短暫離心功能。
- 電動安全門扣、運轉中自動關閉。
- 運轉停止時，自動蜂鳴告示。
- ABS塑鋼外殼並附觀測視窗。
- 自動氣流冷卻循環。
- LED speed and time display.
- The transient centrifugal functions.
- Auto cover lock, auto power off when open.
- Auto buzzer when finish.
- ABS plastics outer covering ; see-through window.
- Auto-cooling air circulation.



Model	Rotor	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)
MCD-2000	REA-1518	14,000	16,000	1.5 / 2.0 x 18	99 : 59	8

Power: AC 110 / 220V , 50 / 60Hz

Size: W290 x D320 x H245 (mm)