

毛細管用離心機 / HEMATOCRIT TYPE CENTRIFUGE


■ HCD-2000 / H-240



HCD-2000



H-240



讀板(Reader)



毛細管(HT tubes)



- 數位式轉速及時間顯示。
- 短暫離心功能。
- 電動安全門扣、運轉中自動關閉。
- 運轉停止時，自動蜂鳴告示。
- ABS塑鋼外殼並附觀測視窗。
- 自動氣流冷卻循環。
- 手動式煞車裝置可迅速減速。
- 回歸式定時設定裝置，時間到自動停止並有鈴聲告示。
- 電源指示燈及保險絲安全裝置。
- 自動氣流冷卻循環。
- 可供觀察校正用之觀測視窗。
- 低噪音結構設計，運轉安靜。

- LED speed and time display.
- The transient centrifugal functions.
- Auto cover lock, auto power off when open.
- Auto buzzer when finish.
- ABS plastics outer covering ; see-through window.
- Auto-cooling air circulation.
- Manual-braking system provides rapid deceleration.
- Regressive 15 min. timer with ring.
- Accessory: Hematocrit reader (order separately)
- Power pilot light and fuse safety protection.
- Auto-cooling air circulation.
- See-through window.
- Low noise.

Model	Rotor	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)
HCD-2000	RHT-240	12,000	15,000	Hematocrit x 24	99:59	10
H-240	RHT-240	Fixed 12,000	15,000	Hematocrit x 24	0 ~ 15 min.	8kg

Power: AC 110 / 220V , 50 / 60Hz

Size: W290 x D320 x H245 (mm)

讀板、毛細管另購。Reader and HT tubes not included.