

油試驗用離心機 / OIL TEST TYPE CENTRIFUGE


■ COD-1004 / 1008

本機適用於測定燃油中水分及沈澱物或潤滑油中溶劑不溶物之用。

- 數位式設定及顯示轉速、離心力、時間。
- 不平衡偵測及自動斷電安全裝置。
- 門蓋安全鎖功能。
- 10組程式記憶功能，可設定不同操作程式。
- 轉速校驗視窗，可供定期校驗轉速。

This machine can be applied to test the water content and sediment in the oil or the insolubility ingredient in the lubricant.

- Microprocessor controlled can be operated precisely.
- LED speed, centrifugal force and time display.
- Auto power off for rotor imbalance.
- Auto cover lock when running.
- 10 memory programs for preset parameters.
- See-through window.



Model	Max. Speed (rpm)	Max. RCF	Oil test bottle (pcs)	Size (mm)
COD-1004	2,000	1,030	4	690 x 770 x 1040
COD-1008		1,120	8	

Power: AC 110 / 220 V , 50 / 60 Hz

試瓶另購。 Tubes not included.